




Throughout human history, the number 3 has always had a unique significance. In Ancient Greek philosophy, the number 3 was considered as the perfect number, the number of harmony, wisdom and understanding.



Correlation coefficients (r) between MC-80 results and manual classification for the five types of leukocytes and the four primary abnormal cells are high for all leukocytes except for reactive lymphocytes.

Comparison of the MC-80 post-classification results with the optical microscope reference method (Carlo Erba™ stain)

	Neutrophils	Lymphocytes	Monocytes	Eosinophils	Basophils	Sensitive granulocytes	Neoplastic cells	Reactive lymphocytes	Total
Pass rate (%)	81.5	87.2	91.8	93.8	98.4	82.2	81.5	75.1	91.3
r	0.949	0.845	0.866	0.949	0.845	0.857	0.907	0.427	0.689
r ²	0.899	0.714	0.754	0.900	0.748	0.734	0.822	0.190	0.790
Notes	0.851	0.850	0.882	0.826	0.822	0.855	1.256	0.295	0.875
Intercept	0.418	0.153	0.578	0.362	0.168	0.045	0.170	0.806	0.307

Note: See Table 53 for the comparison of MC-80 pre-classification data with reference microscope.
 Abbreviation: NRBCs, nucleated red blood cells.
 *Neutrophils include segmented and band neutrophils.
 **sensitive granulocytes=basophils+eosinophils+neutrophils+reactive lymphocytes.
 ***Neoplastic cells include cells classified by the system as lymphoblasts, myeloblasts, abnormal (lymphomatous) lymphocytes, abnormal premyelocyte, and myeloblasts.

The results between pre-classification and the post-classification MC-80 results are compared too. It's reported, differences between the major normal leukocyte classes were minimal, while the most relevant adjustments due to the

Conclusion
(excerpted from "Abstract - Conclusion")

"Our study highlights the outstanding diagnostic performance of this artificial intelligence-based blood film analyzer for hematology patients with circulating abnormal cells. We appreciated the morphological harmonization of cells observed on the screen and those seen in the microscope."



Clinical sensitivity: Identification of films with pathological cells

Stain: May-Grünwald-Giemsa (RAL)



A total of 445 samples were used, and WBC pre-classification values with the MC-80 and DM9600 were compared with the manual reference (microscope), Mindray BC-6800Plus differentials (only normal samples), and confirmed or reclassified images (post-classification).³

Comparative study and statistical analysis

The pre-classification and post-classification of WBC using the MC-80 and DM9600 in normal samples were compared with the manual reference method using the microscope and with the differential leukocyte count (DLC) obtained on the BC-6800Plus.

TABLE 3 Sensitivity, specificity, and predictive values of the MC-80 in identifying films with specific types of cells that do not normally circulate in the blood of healthy subjects.

Active lymphocytes (>1/100 WBC)	NRBCs (>1/100 WBC)		Granulocyte precursors (>1%)	Neoplastic cells* (>0%)	Ri ly (>1%)
3.6%	97.5%	Sensitivity	98.8%	83.8%	97.8%
7.2%	98.9%	Specificity	97.7%	96.1%	97.8%
11.1%	81.4%	PV of positives	83.5%	80.6%	81.1%
17.8%	96.9%	PV of negatives	99.2%	88.4%	97.8%
15.2%	92.0%	Total efficiency	92.0%	85.2%	92.0%

*Neoplastic cells include cells classified by the system as lymphoblasts, atypical lymphocytes, abnormal lymphocytes, lymphoblasts, and atypical lymphocytes. Abbreviations: NRBCs, nucleated red blood cells; PV, predictive value.



[4]Khongjaroensakun N, Chaothai N, Chamchomdao L, Suriyachand K, Paisooksantivatana K. White blood cell differentials performance of a new automated digital cell morphology analyzer: Mindray MC-80. *Int J Lab Hematol*. 2023 Oct;45(5):691-699. doi: 10.1111/ijlh.14119. Epub 2023 Jun 20. PMID: 37338111.

The number 3 holds great inspiration. 3 esteemed morphologists, through their remarkable work on MC-80's 3-year evaluation, have instilled incredible confidence in the analyzer. This has established MC-80 as a trusted tool for serving laboratories worldwide.

References

[1]<https://www.htworld.co.uk/insight/opinion/digital-morphology-what-advances-mean-for-modern-laboratories-and-why-continued-innovation-is-imperative-hm23/>

[2]Zini G, Mancini F, Rossi E, Landucci S, d'Onofrio G. Artificial intelligence and the blood film: Performance of the MC-80 digital morphology analyzer in samples with neoplastic and reactive cell types. *Int J Lab Hematol*. 2023 Dec;45(6):881-889. doi: 10.1111/ijlh.14160. Epub 2023 Aug 28. PMID: 37641457.

[3]Merino A, Laguna J, Rodríguez-García M, Julian J, Casanova A, Molina A. Performance of the new MC-80 automated digital cell morphology analyser in detection of normal and abnormal blood cells: Comparison with the CellaVision DM9600. *Int J Lab Hematol*. 2023 Sep 25. doi: 10.1111/ijlh.14178. Epub ahead of print. PMID: 37746889.



Follow Mindray on Social Media



www.mindray.com

P/N: ENG-Quick Facts About Erythrocyte Sedimentation Rate (ESR) Part A-210X4P-20231010
©2023 Shenzhen Mindray Bio-Medical Electronics Co.,Ltd. All rights reserved.