

# Operating Instructions

## Datascope Duo™





---

*Duo™ is a U.S. trademark of Mindray DS USA, Inc.*

*Navigator™ is a U.S. trademark of Mindray DS USA, Inc.*

*Masimo SET®, LNOP® and CleanShield® are U.S. registered trademarks of Masimo Corp.*

*Nellcor® and OxiMax® are U.S. registered trademarks of Nellcor Puritan Bennett Inc.*

*Copyright © Mindray DS USA, Inc., 2008. All rights reserved. Contents of this publication may not be reproduced in any form without permission of Mindray DS USA, Inc.*



Introduction	-
Getting Started	-
Using the Device	-
Charging the Device	-
Using the Device with a Computer	-
Using the Device with a Smartphone	-
Using the Device with a Tablet	-
Using the Device with a Smart TV	-
Using the Device with a Smart Home System	-
Using the Device with a Smart Car	-
Using the Device with a Smart Watch	-
Using the Device with a Smart Home Hub	-
Using the Device with a Smart Home Camera	-
Using the Device with a Smart Home Lock	-
Using the Device with a Smart Home Light	-
Using the Device with a Smart Home Plug	-
Using the Device with a Smart Home Switch	-
Using the Device with a Smart Home Thermostat	-
Using the Device with a Smart Home Doorbell	-
Using the Device with a Smart Home Door Lock	-
Using the Device with a Smart Home Doorbell Camera	-
Using the Device with a Smart Home Doorbell Camera Viewfinder	-
Using the Device with a Smart Home Doorbell Camera Viewfinder (Optional)	-



**Warnings**

- WARNING:** The Duo monitor is not intended for unsupervised, continuous monitoring. It is for spot-check use only.
- WARNING:** Maintain extreme caution when a defibrillator is in use, avoiding contact with any part of the patient, table or monitor.
- WARNING:** Route cables neatly. Ensure cables, hoses, and wires are away from patient's neck to avoid strangulation. Keep floors and walkways free of cables to reduce the risk of tripping.
- WARNING:** This monitor is not intended for use in an MR environment.
- WARNING:** The Duo monitor is intended for hospital use under the direct supervision of a licensed health care practitioner.
- WARNING:** Do not clean the monitor while it is ON and/or connected to AC power.
- WARNING:** The Duo should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the Duo should be observed to verify normal operation in the configuration in which it will be used.
- WARNING:** Operation of the Duo below the minimum amplitude or value of patient physiological signal may cause inaccurate results.
- WARNING:** Use of accessories, transducers, and cables other than those specified in the manual may result in increased Electromagnetic Emissions or decreased Electromagnetic Immunity of the Duo. It can also cause delayed recovery after the discharge of a cardiac defibrillator.
- WARNING:** Do not use a damaged or broken unit or accessory. Periodically, check all cables (e.g., AC line cord and patient connection cables) for damage that may occur through normal use. Replace cable if damaged in any way.



- CAUTION:** Do not place the SpO<sub>2</sub> sensor on an extremity with an invasive catheter or blood pressure cuff in place.
- CAUTION:** The use of portable and mobile RF communications equipment, in the proximity of the Duo, can affect the performance of this monitor.
- CAUTION:** Use only Mindray DS accessories with this product. For a comprehensive listing of Duo Accessories refer to section 4.0, "Accessories."
- CAUTION:** The patient size selection should be matched to the actual patient before monitoring begins.
- CAUTION:** Tissue damage or inaccurate measurement may be caused by incorrect SpO<sub>2</sub> sensor application or use, such as wrapping too tightly, applying supplemental tape, failing to inspect the sensor site periodically or failing to position appropriately. Carefully read the SpO<sub>2</sub> sensor directions and all precautionary information before use.
- CAUTION:** Excessive ambient light may cause inaccurate SpO<sub>2</sub> measurements. In such cases, cover the sensor site with opaque material.
- CAUTION:** The cuff must be properly applied to the patient's limb before inflating. If it is inflated without being securely wrapped, damage to the cuff can result.
- CAUTION:** This product contains natural rubber latex which may cause allergic reactions. This refers specifically to the large adult gray blood pressure cuff (0998-00-0003-35).
- CAUTION:** If the device is accidentally saturated with any liquid, immediately discontinue use and contact service personnel.

- CAUTION:** Inaccurate SpO<sub>2</sub> measurements may be caused by:
- incorrect sensor application or use
  - significant levels of dysfunctional hemoglobins, (e.g., carboxyhemoglobin or methemoglobin)
  - intra-vascular dyes such as indocyanine green or methylene blue
  - exposure to excessive illumination such as surgical lamps (especially ones with a xenon light source), bilirubin lamps, fluorescent lights, infrared heating lamps, or excessive ambient light. In such cases, cover the sensor site with opaque material.
  - excessive patient movement
  - venous pulsations
  - electro-surgical interference
  - placement of a sensor on an extremity that has a blood pressure cuff, arterial catheter or intra-vascular line.
  - nail polish or fungus
- CAUTION:** In certain situations in which perfusion and signal strength are low, such as in patients with thick or pigmented skin, inaccurately low SpO<sub>2</sub> readings will result. Verification of oxygenation should be made, especially in patients with chronic lung disease, before instituting any therapy or intervention.
- CAUTION:** Many patients suffer from poor peripheral perfusion due to hypothermia, hypovolemia, severe vasoconstriction, reduced cardiac output, etc. These symptoms may cause a loss in vital sign readings.
- CAUTION:** If the SpO<sub>2</sub> sensor or patient cable are damaged in any way, discontinue use immediately. To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize.
- CAUTION:** When applying the SpO<sub>2</sub> sensor to the patient, ensure proper positioning, alignment and skin integrity. Exercise extreme caution with poorly perfused patients.
- CAUTION:** When equipped with Masimo SpO<sub>2</sub>, use only Masimo oxygen sensors and cables. Use of other oxygen sensors may cause improper oximeter performance.
- CAUTION:** When equipped with Nellcor SpO<sub>2</sub>, use only Nellcor oxygen sensors and cables. Use of other oxygen sensors may cause improper oximeter performance.
- CAUTION:** Use only Mindray DS blood pressure cuffs and hoses with the Duo.

- 
- CAUTION:** A patient's skin is sometimes fragile (i.e., on pediatric and geriatric patients, or due to physiological conditions). In these cases, a longer duration between NIBP measurements should be considered to decrease the number of cuff inflations over a period of time. In extreme cases, a thin layer of soft roll or cotton padding may be applied to the limb in order to cushion the skin when the cuff is inflated. This may affect NIBP performance and should be used with caution.
- CAUTION:** Please consult a physician for interpretation of blood pressure measurements.
- CAUTION:** A blood pressure measurement can be affected by the position of the patient, and his/her physiological condition as well as other factors, such as patient movement.
- CAUTION:** Any condition that may affect the regularity and strength of arterial pressures (such as patient movement, cardiac arrhythmias, restriction of hose, etc.), will affect the accuracy and ability to measure the NIBP.
- CAUTION:** When cleaning SpO<sub>2</sub> sensors, do not use an excessive amount of liquid. Wipe the sensor surface with a soft cloth, dampened with a cleaning solution.
- CAUTION:** Do not subject the SpO<sub>2</sub> sensor to autoclaving.
- CAUTION:** Do not use SpO<sub>2</sub> sensors or cables that are damaged or have deteriorated.
- CAUTION:** Some disinfectants may cause skin irritation. Please rinse the NIBP cuffs thoroughly with water to remove any residual disinfectants.
- CAUTION:** Using dark colored soaps may stain the NIBP cuffs. Test a single cuff to ensure that no damage will occur.
- CAUTION:** Disposable NIBP cuffs can be cleaned using a mild soap solution and dried with a clean cloth.
- CAUTION:** Replace the Lithium Ion battery with part number 0146-00-0079 only.
- CAUTION:** Remove the battery if the Duo is not likely to be used for an extended period of time.
- CAUTION:** Remove the battery prior to shipping the Duo.

7-11

- NOTE:** Potential hazards due to errors in software or hardware have been minimized by actions taken in accordance with IEC 60601-1-4.
- NOTE:** Information codes and error codes with corresponding explanations are provided to assist in the identification and correction of problems that may occur with the monitor.
- NOTE:** The comparison testing conducted via the auscultatory method used both Phase 4 and Phase 5 Korotkoff sounds. A report of the study finding for the auscultatory method is available by contacting Technical Support (201) 995-8116.
- NOTE:** The use of this equipment is restricted to one patient at a time.



















Figure 1-1


SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Warning: Read the instructions carefully before using the device.		Person: This symbol indicates a user action or a person-related warning.
	High Voltage: Caution against touching live parts.		Person with High Voltage: Warning of high voltage danger to the person.
	Grounding: The device must be grounded.		Medical: Caution against using the device on children or pregnant women.
	Fire: The device is flammable.		No Fire: Do not use open flames.
	No Open Flame: Do not use open flames.		No Person: Do not use the device on children or pregnant women.
	No Person: Do not use the device on children or pregnant women.		No Person with High Voltage: Do not use the device on children or pregnant women.
	No Person with High Voltage: Do not use the device on children or pregnant women.		No Person with High Voltage: Do not use the device on children or pregnant women.
<p>For more information, see the Duo™ Operating Instructions (11) 1.</p>			

**This page intentionally left blank.**

---


1.0 *General Product Description*

---

1.1 

## 1.2

### 1.2.1




**FIGURE 1-1**

#### **1. Digital Display**

**Parameter Tiles**

The parameter tiles display vital signs and other patient data. The tiles are arranged in a vertical column on the right side of the screen. Each tile contains a numerical value and a unit of measurement. The tiles are color-coded: red for NIBP, green for Pulse Rate, and blue for SpO<sub>2</sub>. The tiles are labeled with letters a, b, and c.



**FIGURE 1-2** Parameter Tiles

**a. NIBP**

The NIBP tile displays the following information:

- Sys.** (Systolic Blood Pressure)
- Dia.** (Diastolic Blood Pressure)
- MAP** (Mean Arterial Pressure)

Units: mmHg, kPa.

**b. Pulse Rate**

The Pulse Rate tile displays the following information:

- Pulse rate (bpm)
- Heart rate icon
- Source icon

Units: bpm.

**c. SpO<sub>2</sub> (Non-invasive)**

The SpO<sub>2</sub> tile displays the following information:

- SpO<sub>2</sub> value (%)
- SpO<sub>2</sub> icon

Units: %.

**d. Information Codes**


The information codes are displayed in the bottom right corner of the screen. They provide additional details about the patient's status and the device's settings.

**LED Indicators**




**FIGURE 1-3**

**a. Battery Charging**




**b. AC Power**





**c. Patient Size**



**d. NIBP Start/Stop**




## 2. Keypad




**FIGURE 1-4**


### a. Power ON/OFF




**NOTE:** If the Power ON/OFF key is depressed for less than two (2) seconds, the monitor will not power OFF.




### b. Clear/Next Patient




### c. Patient Size



### d. NIBP Start/Stop




### 3. Connector Panel




**FIGURE 1-5**

**a. NIBP Pneumatic Fitting**





**b. SpO<sub>2</sub> Receptacle (optional)**




\* *Quick Connect Pneumatic Fittings available from Rectus-TEMA Corporation.*

1.2.2




**FIGURE 1-6**

**1. Serial Port**




**2. Equipotential Lug**




**NOTE:** Ensure that when connecting external devices to the unit all equipotential terminals are connected.

**3. AC Receptacle**



**NOTE:** The power supply, and battery charger (if the battery is installed) are active any time AC power is supplied, regardless of whether the monitor is ON or OFF.

**4. Mounting Alignment Slot**







When the Duo is in Standby Mode, it will automatically return to Normal Monitoring Mode if the battery voltage drops below a certain level. This is to ensure that the Duo is always ready to monitor the system.

EVENT	TIME PERIOD TO ENTER STANDBY MODE
Power ON/OFF	1.0 hr.
Power ON/OFF	1.0 hr.
Power ON/OFF	1.0 hr.
Power ON/OFF	1.0 hr.

**Power ON/OFF.** The Duo will automatically return to Normal Monitoring Mode if the battery voltage drops below a certain level. This is to ensure that the Duo is always ready to monitor the system.

**NOTE:** The Duo cannot be powered OFF while in Standby Mode. It must return to Normal Monitoring Mode before it can be powered OFF.

### 2.1.3

The Duo will automatically return to Normal Monitoring Mode if the battery voltage drops below a certain level. This is to ensure that the Duo is always ready to monitor the system.

### 2.1.4

The Duo will automatically return to Normal Monitoring Mode if the battery voltage drops below a certain level. This is to ensure that the Duo is always ready to monitor the system.

## 2.2


1. Turn the Duo™ on. (See Figure 2-1.)

1. Turn the Duo™ on. (See Figure 2-1.)
2. Press the power button (see Figure 2-2).



**FIGURE 2-1** Turn the Duo™ on.


2. Press the power button (see Figure 2-2).




**FIGURE 2-2** Press the power button.

3. Press the power button again to turn the Duo™ off. (See Figure 2-3.)


**NOTE:** The Lithium Ion battery is shipped in a partially charged state and must be fully charged prior to its first use.



3. 



**NOTE:** Optimum battery runtime is achieved after 3 charge/discharge cycles.

4. 

### 2.2.1



Duo  N<sub>1</sub> (%) 

 mmHg  mmHg.



## 2.3

... Duo.

### 2.3.1

...

**CAUTION:** A patient's skin is sometimes fragile (i.e., on pediatric and geriatric patients, or due to physiological conditions). In these cases, a longer duration between NIBP measurements should be considered to decrease the number of cuff inflations over a period of time. In extreme cases, a thin layer of soft roll or cotton padding may be applied to the limb in order to cushion the skin when the cuff is inflated. This may affect NIBP performance and should be used with caution.

**CAUTION:** Please consult a physician for interpretation of blood pressure measurements.


**CAUTION:** A blood pressure measurement can be affected by the position of the patient, and his/her physiological condition as well as other factors, such as patient movement.




**CAUTION:** Any condition that may affect the regularity and strength of arterial pressures (such as patient movement, cardiac arrhythmias, restriction of hose, etc.), will affect the accuracy and ability to measure the NIBP.

... Duo t. ... NIBP Start/Stop ... Duo ...

PATIENT SIZE SETTING	DEFAULT CUFF INFLATION PRESSURE
...	...
...	...


... Duo ... MAP ...






**NOTE:** Pressing the NIBP Start/Stop key while the NIBP measurement is in progress will stop the measurement and deflate the cuff.


### 2.3.2



**NOTE:** If the optional SpO<sub>2</sub> is not purchased with the Duo, the Pulse Rate source will be NIBP by default.




**NOTE:** Pulse Rate from NIBP is a static value since NIBP is a static, one-time measurement.










### 2.3.3.1




**CAUTION:** When equipped with Masimo SpO<sub>2</sub>, use only Masimo oxygen sensors and cables. Use of other oxygen sensors may cause improper oximeter performance.


**NOTE:** Refer to instructions included with each SpO<sub>2</sub> sensor and cable for proper placement and use.

1. 
2. 
3. 

### 2.3.3.2





**CAUTION:** When equipped with Nellcor SpO<sub>2</sub>, use only Nellcor oxygen sensors and cables. Use of other oxygen sensors may cause improper oximeter performance.




**NOTE:** Refer to instructions included with each SpO<sub>2</sub> sensor and cable for proper placement and use.

1. Plug the SpO<sub>2</sub> cable into the Duo.
2. Turn the Duo on.
3. Place the SpO<sub>2</sub> sensor on the patient's finger.
4. The Duo will display the SpO<sub>2</sub> reading.



**NOTE:** To disconnect the cable from the Duo, squeeze the tabs on the sides of the connector and then pull it straight out.

## 2.4



**NOTE:** Information codes and Error codes that are marked with an asterisk (\*) are one-time codes that can be cleared from the display by pressing the Clear/Next Patient key.

2.4.1

MESSAGE TYPE	CODE	DESCRIPTION	REASON
M	*NN		
	*		
	*		
N	1		
	10		
	11		
	1		
	1		
	1		

MESSAGE TYPE	CODE	DESCRIPTION	REASON
M << N, P N1	P N1		
	P N1		
	P N1	T R	
	P N1	P F M <	
	P N1	N N	
	P N1	N T	
F M < N, P N1		N F N	


## 2.4.2

1 1 1 1 1 1


MESSAGE TYPE	CODE	DESCRIPTION	REASON
M	1	M P < P . "	M P < P . "
	1	M P . N "	M P . N "
		<	<

---


3.1



## 3.2




**WARNING:** Do not clean the monitor while it is ON and/or connected to AC power.



## 3.3



### 3.3.1



**NOTE:** Refer to the individual instruction sheets that are packaged with each sensor.

1. 1.
2. 2.

**CAUTION:** When cleaning SpO2 sensors, do not use an excessive amount of liquid. Wipe the sensor surface with a soft cloth, dampened with a cleaning solution.

**CAUTION:** Do not subject the SpO2 sensor to autoclaving.

**CAUTION:** If the SpO2 sensor or patient cable are damaged in any way, discontinue use immediately. To prevent damage, do not soak or immerse the sensor in any liquid solution. Do not attempt to sterilize.

**CAUTION:** Do not use SpO2 sensors or cables that are damaged or have deteriorated.

### 3.3.2

**NOTE:** Accuracy of cuff-pressure transducers/indicators is to be verified at intervals specified by the manufacturer.

#### 3.3.2.1


##### Cleaning

**NOTE:** Machine washing may shorten the service life of the cuff.

##### Disinfection


**NOTE:** Prolonged use of disinfectant may cause discoloration of the cuff.

1. Wash the cuff with fresh water.
2. Rinse the cuff with fresh water.
3. Disinfect the cuff with 70% isopropanol.
4. Dry the cuff with a clean, dry cloth.




**CAUTION:** Do not dry clean the cuff.  
Do not press the cuff with a hot iron.  
Do not use detergent and disinfectant other than fresh water or 70% isopropanol.  
Clean and disinfect the cuff according to the instructions.

### 3.3.2.2




### 3.3.2.3




**NOTE:** Disposable cuffs can be cleaned using a mild soap solution and dried with a clean cloth. For Cuffs with bladders, remove bladder before cleaning.

### 3.4




**CAUTION:** Replace the Lithium Ion battery with part number 0146-00-0079 only.

- 1.
- 2.





**FIGURE 3-1**

- 3.




**FIGURE 3-2**

4. 
5. 



**CAUTION:** Remove the battery if the Duo is not likely to be used for an extended period of time.

**CAUTION:** Remove the battery prior to shipping the Duo.

• **Duo** 

## Accessories

---


### 4.1

#### Accessories


Accessories are available at <http://www.hewlett-packard.com>, 0020-00-0123-01

DESCRIPTION	PART NUMBERS
(1) <i>[Illegible]</i>	- - - 1
(1) <i>[Illegible]</i>	- - - 1
(1) <i>[Illegible]</i>	(1) - ( - 1


4.2




4.2.1



DESCRIPTION	PART NUMBERS
	- -
	- -



DESCRIPTION	PART NUMBERS
	( - ( - (
	( - 1- (
	( - - (
	( - - (
	( - - (
	( - - (
	( - - (



DESCRIPTION	PART NUMBERS
	( - ( - (
	( - 1- (
	( - - (
	( - - (
	( - - (
	( - - (
	( - - (



**This page intentionally left blank.**









5.2

MAP (mmHg) 100 110 120 130 140 150 160 170 180 190 200


5.2.1



5.2.2




5.2.3



**NOTE:** Mean Arterial Pressure (MAP) is defined as:  
 Mean Pressure 1 = Mean Pressure determined from the oscillometric profile  
 Mean Pressure 2 = (2\*diastolic + systolic) / 3  
 Mean Pressure Displayed = (Mean Pressure 1 + Mean Pressure 2) / 2

5.2.4



5.2.5

Flow Rate (L/min)

Flow Rate (L/min) 0.1-12  
 Flow Rate (L/min) 0.1-12

	ADULT MODE	PEDIATRIC MODE
Flow Rate (L/min)	0.1-12	0.1-12

5.2.6

Flow Rate (L/min) 0.1-12

5.2.6.1

Flow Rate (L/min) 0.1-12

Flow Rate (L/min) 0.1-12

ADULT MODE	PEDIATRIC MODE
Flow Rate (L/min) 0.1-12	Flow Rate (L/min) 0.1-12

Flow Rate (L/min) 0.1-12

ADULT MODE	PEDIATRIC MODE
Flow Rate (L/min) 0.1-12	Flow Rate (L/min) 0.1-12

5.2.6.2

Flow Rate (L/min)

Flow Rate (L/min) 0.1-12

Flow Rate (L/min) 0.1-12

5.2.6.3

Flow Rate (L/min)

Flow Rate (L/min) 0.1-12

5.2.6.4

Flow Rate (L/min)

Flow Rate (L/min) 0.1-12


5.2.6.5

Flow Rate (L/min)

Flow Rate (L/min) 0.1-12


Flow Rate (L/min) 0.1-12

5.2.6.6




PATIENT SIZE	PRESSURE INCREMENT (DEPENDING ON ALGORITHM)	DEFAULT START PRESSURE
1	2	10
2	2	10


5.2.6.7



5.2.7



5.2.7.1



5.2.7.2




Figure 1:

Figure 1: Accuracy of SpO<sub>2</sub> readings

PATIENT SIZE	SATURATION RANGE	
	70% to 100%	0 - 69%
Adult	±1.0%	±0.5%
Child	±1.0%	±0.5%

Figure 2: Accuracy of SpO<sub>2</sub> readings

PATIENT SIZE	SATURATION RANGE	
	70% to 100%	0 - 69%
Adult	±1.0%	±0.5%
Child	±1.0%	±0.5%

Figure 3: Accuracy of SpO<sub>2</sub> readings

PATIENT SIZE	SATURATION RANGE	
	70% to 100%	0 - 69%
Adult	±1.0%	±0.5%
Child	±1.0%	±0.5%

Figure 4: Accuracy of SpO<sub>2</sub> readings

**LOW PERFUSION CONDITIONS**

PULSE AMPLITUDE	% TRANSMISSION	SATURATION ACCURACY	PULSE RATE ACCURACY
≥1%	≥%	±1.0%	±0.5%

Figure 5: Accuracy of SpO<sub>2</sub> readings

Figure 6: Accuracy of SpO<sub>2</sub> readings

Figure 7: Accuracy of SpO<sub>2</sub> readings

Figure 8: Accuracy of SpO<sub>2</sub> readings

PATIENT SIZE	PULSE RATE RANGE	ACCURACY	
		NO MOTION CONDITIONS	DURING MOTION CONDITIONS <sup>1</sup>
Adult	1-180	±1.0%	±1.0%

The table contains several rows of data, likely representing different patient parameters and their specifications. Due to the extreme fading, the specific values and units are not legible. The structure appears to be a list of parameters with associated numerical values and units.

**NOTE:** The sensor measurement wavelengths are nominally 660 nm for the red LED and 940 nm for the infrared LED. Maximum optical power output for LED is 4 mW.

### 5.2.7.3

2

2

%

%

SENSOR	ACCURACY
	%
	%
	%
	%

RANGE	ACCURACY

**NOTE:** The sensor measurement wavelengths are nominally 660 nm for the red LED and 890 nm for the infrared LED. Maximum optical power output for LED is 4 mW.



5.2.11  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

5.2.11.1

$$\begin{aligned} \left( \frac{1}{\tau_{12}} \right) &= \frac{1}{\tau_{12}} \\ \left( \frac{1}{\tau_{21}} \right) &= \frac{1}{\tau_{21}} \\ \left( \frac{1}{\tau} \right) &= \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}} \end{aligned}$$

5.2.11.2

$$\begin{aligned} \left( \frac{1}{\tau} \right) &= \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}} \\ \left( \frac{1}{\tau} \right) &= \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}} \end{aligned}$$

5.2.12  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

5.2.13  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**No**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\left( \frac{1}{\tau} \right) = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**No**  $\left( \frac{1}{\tau} \right) = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\left( \frac{1}{\tau} \right) = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**No**  $\left( \frac{1}{\tau} \right) = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$


**Duo**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**Duo**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**No**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

**F**  $\frac{1}{\tau} = \frac{1}{\tau_{12}} + \frac{1}{\tau_{21}}$

Parameter 1	Value 1
Parameter 2	Value 2
Parameter 3	Value 3
Parameter 4	Value 4



**NOTE:** The Duo needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.

**NOTE:** Portable and mobile RF communications equipment can affect the Duo. See TABLE 5-1 through TABLE 5-4.

**TABLE 5-1**

**GUIDANCE AND DECLARATION - ELECTROMAGNETIC EMISSIONS**

EMISSIONS TEST	COMPLIANCE	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
EN 55011 Class B	EN 55011 Class B	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.
EN 55012 Class B	EN 55012 Class B	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.
EN 61010-1	EN 61010-1	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.

**TABLE 5-2**

**GUIDANCE AND DECLARATION - ELECTROMAGNETIC IMMUNITY**

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
EN 61010-1	EN 61010-1	EN 61010-1	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.
EN 61326-2-2	EN 61326-2-2	EN 61326-2-2	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.
EN 61010-1	EN 61010-1	EN 61010-1	The Duo is designed to be used in an electromagnetic environment where the maximum magnetic field strength is 100 A/m (2.5 mT) at 150 kHz to 1 MHz, and 3 V/m at 1.5 MHz to 100 MHz.

**TABLE 5-2**

**GUIDANCE AND DECLARATION - ELECTROMAGNETIC IMMUNITY**



**Duo** is designed to meet the following immunity test levels:

IMMUNITY TEST	IEC 60601 TEST LEVEL	IEC 60601 TEST LEVEL	<b>Duo</b>
Electrostatic Discharge (ESD)	± 8 kV (air discharge)	± 4 kV (contact discharge)	± 8 kV (air discharge)
Electromagnetic Interference (EMI) - Conducted	30 V (0.5 A)	30 V (0.5 A)	30 V (0.5 A)
Electromagnetic Interference (EMI) - Radiated	3 V/m (100 kHz - 1 GHz)	3 V/m (100 kHz - 1 GHz)	3 V/m (100 kHz - 1 GHz)
Power Frequency Magnetic Field (50/60 Hz)	3 A/m (30 A)	3 A/m (30 A)	3 A/m (30 A)
Power Frequency Electric Field (50/60 Hz)	10 kV/m	10 kV/m	10 kV/m
Surge Immunity	1 kV (10 A)	1 kV (10 A)	1 kV (10 A)
Fast Transient Burst (FTB)	2 kV (10 A)	2 kV (10 A)	2 kV (10 A)
Immunity to RF Fields (100 kHz - 1 GHz)	3 V/m	3 V/m	3 V/m
Immunity to RF Fields (1.5 - 100 MHz)	3 V/m	3 V/m	3 V/m
Immunity to RF Fields (100 MHz - 1 GHz)	3 V/m	3 V/m	3 V/m

*U<sub>i</sub>* is the A.C. mains voltage prior to application of the test level.

**TABLE 5-3**

**GUIDANCE AND DECLARATION - ELECTROMAGNETIC IMMUNITY**

IMMUNITY TEST	IEC 60601 TEST LEVEL	COMPLIANCE LEVEL	ELECTROMAGNETIC ENVIRONMENT - GUIDANCE
			<p><b>Recommended separation distance</b></p> $d = 0.1 \sqrt{P}$ $d = 0.1 \sqrt{P} \quad V, F$ $d = 1. \sqrt{P} \quad V, F$
<p><math>P</math> is the power of the transmitter in Watts (W).  <math>d</math> is the recommended separation distance in meters (m).  </p>			

**NOTE:** At 80 MHz and 800 MHz, the higher frequency range applies.

**NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

*a* Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Duo is used exceeds the applicable RF compliance level above, the Duo should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Duo.

*b* Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

**TABLE 5-4**


**RECOMMENDED SEPARATION DISTANCES BETWEEN PORTABLE AND MOBILE RF COMMUNICATIONS EQUIPMENT AND THE DUO**

RATED MAXIMUM OUTPUT POWER OF TRANSMITTER W (WATTS)	SEPARATION DISTANCE ACCORDING TO FREQUENCY OF TRANSMITTER M (METERS)		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$M = 0.1 \sqrt{P}$	$M = 0.1 \sqrt{P}$	$M = 1.0 \sqrt{P}$
0.1	0.1	0.1	0.1
1	0.1	0.1	0.1
10	0.1	0.1	0.1
100	0.1	0.1	0.1
1000	0.1	0.1	0.1

**NOTE:** At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

**NOTE:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.


**NOTE:** The Duo is intended for use in the electromagnetic environment specified below. The customer or the user of the Duo should assure that it is used in such an environment.




**NOTE:** Mains power quality should be that of a typical commercial or hospital environment. If the user of the Duo requires continued operation during power mains interruptions, it is recommended that the Duo be powered from an uninterruptible power supply or a battery.




### 5.3





# 5.4



- a.
- b.
- c.



