

## TM70/TM80/BP10 Default Settings

Hospital, City, ST \_\_\_\_\_ Department \_\_\_\_\_

Clinical/Sales/Service Representatives \_\_\_\_\_

Date \_\_\_\_\_ Software Version \_\_\_\_\_ Signature, Title \_\_\_\_\_

Maintenance				
General	Display Auto Off: 1 min <input type="checkbox"/> 2min <input type="checkbox"/> 5min <input type="checkbox"/> 15 min <input type="checkbox"/> 30min <input type="checkbox"/> Off <input type="checkbox"/>			
	Language: English			
	Device Location:		Facility:	
	Department: Fixed <input type="checkbox"/> Unfixed <input type="checkbox"/>		Department:	
	Patient Location	Bed/Room: Fixed <input type="checkbox"/> Unfixed <input type="checkbox"/>	Room No.:	Bed No.:
	Patient Info Display: Auto <input type="checkbox"/> Patient Name <input type="checkbox"/> Patient ID <input type="checkbox"/> Bed No. <input type="checkbox"/>			
	Notch Filter: 60 <input type="checkbox"/> 50 <input type="checkbox"/> Off			
	ECG lead label : AHA <input type="checkbox"/>			
	SpO2 Massimo <input type="checkbox"/> Nellcor <input type="checkbox"/>			
	SpO2 tone: Mode1 <input type="checkbox"/> Mode 2 <input type="checkbox"/>			
	Date/Time	Time: 12 hr <input type="checkbox"/> 24 hr <input type="checkbox"/>	Date: DD/MM/YYYY <input type="checkbox"/>	YY/MM/DD <input type="checkbox"/> MM/DD/YY <input type="checkbox"/>
	Resp	On <input type="checkbox"/> Off <input type="checkbox"/>		
Alarms	Style ISO <input type="checkbox"/> Mode1 <input type="checkbox"/> Mode2 <input type="checkbox"/>			
	Minimum Alarm Volume			
	Reminder tone On <input type="checkbox"/> Off <input type="checkbox"/>			
	Reminder Interval 1 min <input type="checkbox"/> 2min <input type="checkbox"/> 3min <input type="checkbox"/> 5min <input type="checkbox"/> 10min <input type="checkbox"/>			
	ECG Lead off Low <input type="checkbox"/> Med <input type="checkbox"/> High <input type="checkbox"/>			
	SpO2 sensor off Low <input type="checkbox"/> Med <input type="checkbox"/> High <input type="checkbox"/>			
	Alarm Reset By Other Bed: On <input type="checkbox"/> Off <input type="checkbox"/>			
Quick Keys				
Screen lock	Off <input type="checkbox"/> Locked <input type="checkbox"/> View Only <input type="checkbox"/>			
Auto Enter Mon. Mode	No Central Monitoring <input type="checkbox"/> Lethal Arrhy. Alarm <input type="checkbox"/>			
Alarms				
<b>*See alarm report printout</b>				
Review				
<b>Tabular Trends</b>	<b>Interval:</b>			
Display Setup				
Portrait		Landscape		
Rows 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/>		Rows 2 <input type="checkbox"/> 3 <input type="checkbox"/> 4 <input type="checkbox"/>		
Order		Order		
Default Orientation	Portrait <input type="checkbox"/> Landscape <input type="checkbox"/>			
Brightness				
Audio volume	Tech alarm:	Touch screen click:	Systole Beep:	



<b>TD60 Default Settings</b>		<b>0002-08-37031</b>
		<b>1.0</b>
A	CR 42661 01/29/2016	Original Release as Rev A
B	CR 44309 07/06/2018	Updated to Rev B
1.0	CR2309-00291 10/05/2023	Updated to Rev 1.0. Text changes throughout the PowerPoint.